

## **Sizing the TAFO:**

1. **Fit it. Measure for an accurate size.** Allow for a hand's width between back of knee and top of splint.\*

\*If fitting the TAFO *to attach to the PKO*, measure from the base of the heel to the crease behind the knee (popliteal space) with the leg extended.

2. Toes should not extend over the splint.
3. The splint must allow full knee flexion.
4. Do not over tighten the straps. There should be space for a finger to slide between the skin or dressing and the strap pad.

## **TAFO Item #'s, Sizes and Specifications:**

**100/ 2100 Xlg: 20-23 inches (50.8-58.4 cm) calf length.**

Fits a tall person ~6'1 1/2" (182.5 cm) or taller and accommodates for feet size ≥13

**101/ 2101 Lg: 18-20 inches (45.7-50.8 cm) calf length.**

Fits a person ~5'10"- 6'1" (173 - 182 cm) tall with average size feet (10-12.5)

**102/ 2102 Med: 16-18 inches (40.6-45.7 cm) calf length.**

Fits a person ~5'5"- 5'9 1/2" (162.5 - 172.5 cm) tall with average size feet (7-9.5)

**103/ 2103 Sm: 14-16 inches (35.6-40.6 cm) calf length.**

Fits a person ~5'0"- 5'4 1/2" (150 - 160 cm) tall with average size feet (5-6.5 or 7M)

**104 Youth: 12 1/2-14 inches (31.8-35.6 cm) calf length.**

Fits a small person who is ~3'5" - 5'0" (102.5 - 149 cm) tall with feet size (3-5)

**105 Child: 9-12 inches (22.9-30.5 cm) calf length.**

**106 Toddler: 5-8 inches (12.7-20.3 cm) calf length.**

Disclaimer: PFS Med orthoses are designed to be sized and customized by a qualified practitioner. PFS Med, Inc does not express or imply responsibility for injury or impairment in skin integrity resulting from splints improperly fit or inadequately supervised.